



Product Data Sheet

Medical Tap Filter Fixed Position 31 Day

t-safe.com

 **T safe**

Medical Tap Filter Fixed Position 31 Day



Features

- Sterilising grade 0.2 µm rated membrane filter
- Up to 31 day life cycle
- Compatible with swan neck and wall mounted taps
- Laminar Flow outlet reduces splashing
- Constant flow technology ensures optimal flow performance over life cycle
- Antimicrobial additive protection reduces risk of retrograde contamination
- Integral check valve prevents backflow contamination upstream
- Traceability with GS1 Datamatrix, enabling electronic tracking of filters

Order Information

Item	Code
Medical Tap Filter Fixed Position 31 Day	02-831401
Medical Tap Filter Fixed Position 31 Day (sterile packed)	02-831402

Technical Data

Max. operating pressure	5 bar
Max. operating temperature	60°C
Flow rate	6 L/min
Max. filter life cycle	31 days*
Connection	Adapter
Colour	White
Microbial retention	
Sterilising grade membrane filter, effective against all waterborne pathogens. Independently tested and validated in accordance with ASTM F838 against:	
<ul style="list-style-type: none">• <i>Legionella pneumophila</i>• <i>Pseudomonas aeruginosa</i>• <i>Brevundimonas diminuta</i>	
Microbial reduction	
T-safe antimicrobial additive is proven to reduce surface bacteria contamination by 99.9%. Independently tested in accordance with ISO 22196.	

Approvals

ASTM F838, ISO 22196, WRAS, DVGW W270. Compliant with EU food and drinking water directives.

* Dependent on local water quality and usage.

All information contained within this data sheet is based on the current knowledge and experience of T-safe. It is provided in good faith and is intended to serve only as a guideline for the selection and use of products. T-safe products should always be handled in accordance with the applicable Instruction for Use advice. The current data sheet is available at www.t-safe.com.